



# Escape-Mattress Bariatric



## For the evacuation of obese individuals

### What's an Escape-Mattress?

The **Escape-Mattress Bariatric** is designed for evacuating obese people (up to 300kg) in lying position. The **Escape-Mattress Bariatric** needs to be operated by several trained persons. Also suitable for spiral staircases.

The **Escape-Mattress Bariatric** can be used to move evacuees along the floor (by using at least 4 persons and using the additional towing straps) and have 12 handles placed on the sides to allow lifting and easy maneuverability.

### Advantages of the Escape-Mattress Bariatric:

- TÜV certified on 300KG.
- Extra wide version.
- Made of high quality materials (fire resistance classification M2).
- Extra-long towing straps to help operators with improved handling practices.
- Handgrips to assist with sideways manoeuvrability.
- Carrying people is no longer necessary.
- Ready for use within seconds.
- Safe evacuation down the stairs.
- Speed of descent is determined by the user.
- Fits into your emergency evacuation plan.
- CE, TÜV NORD and ISO tested.

### Specifications Escape-Mattress Bariatric

Article number:	21080
Dimensions (length x width x thick):	± 200 x 86 x 2 cm
Maximum payload:	300 kg
Weight:	±10 kg incl. bag
Equipped with additional, movable drag straps:	✓
Ready for use within seconds:	✓
Speed of descent to be determined by the user:	✓
Blood, urine & acid resistant:	✓
Heat and flame retardant:	✓
Fire resistance classification M2:	✓
Equipped with headrest, body-covers & foot sack:	✓
Including storage bag:	✓
User manual:	✓
Photoluminescent identity sign:	✓
CE, TÜV NORD certified & ISO tested:	✓
Warranty:	2 years

### Our advice:

Make sure you have an **Escape-Mattress Bariatric** available, so you can help obese people who rely on you to be rescued safely. Even when you need to tackle spiral staircases.

**"Failure to prepare for an emergency has its consequences."**



**E:** [info@escape-mobility.com](mailto:info@escape-mobility.com)

**Escape Mobility Company Ltd.**  
85 Great Portland Street, 1st floor,  
London, W1W 7LT

**T:** +44 (0)1908 229 028

**W:** [www.escape-mobility.com](http://www.escape-mobility.com)

**Escape Mobility International BV.**  
Jan Campertstraat 7 – Unit 1.14,  
6416 SG, Heerlen.

**T:** +31 (0)45 528 0670

...your way to safety

